

Mouse anti-Myeloid Specific Antigen, clone BM-5 (Monoclonal)

Clone no. BM-5

MONOSAN

---

Product name	Mouse anti-Myeloid Specific Antigen, clone BM-5 (Monoclonal)
Host	Mouse
Applications	IHC-fr
Species reactivity	human
Conjugate	-
Immunogen	nuclei from Pokeweed mitrogen stimulated human peripheral blood lymphocytes
Isotype	IgG1
Clonality	Monoclonal
Clone number	BM-5
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Myeloid Specific Antigen, clone BM-5 (Monoclonal)

Clone no. BM-5

MONOSAN

**Additional info**

BM-5 is specific marker for human myeloid cells and an early marker of myeloid differentiation. It recognizes a nuclear and cytoplasmic antigen present in granulocytes, monocytes, and myeloid precursor cells. It also reacts with a subset of myeloid leukemia cells. BM-5 has no reactivity with any other cell type in human tissues.

**References**

1. Epstein, AL et al. Blood 70:1124-1130 (1987)
2. Murao et al, Cancer Research 45: 791-795 (1985)
3. -
4. -
5. -

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**