

Mouse anti-Phosphotyrosine, clone P9V6 (Monoclonal)

Clone no. P9V6

MONOSAN

Product name	Mouse anti-Phosphotyrosine, clone P9V6 (Monoclonal)
Host	Mouse
Applications	ELISA,WB
Species reactivity	n/a
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	P9V6
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Phosphotyrosine, clone P9V6 (Monoclonal)

Clone no. P9V6

MONOSAN

Additional info

This highly specific monoclonal antibody P9V6 reacts with phosphotyrosine. It does not cross react with phosphoserine or phosphothreonine.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES