

Mouse anti-Phosphotyrosine, clone P-TYR-01 (Monoclonal)

Clone no. P-TYR-01

MONOSAN

---

Product name	Mouse anti-Phosphotyrosine, clone P-TYR-01 (Monoclonal)
Host	Mouse
Applications	WB, ICC, FC
Species reactivity	Species independent
Conjugate	-
Immunogen	Phosphotyrosine conjugated to bovine serum albumin.
Isotype	IgG1
Clonality	Monoclonal
Clone number	P-TYR-01
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON2055

**MONOSAN**<sup>®</sup>

Mouse anti-Phosphotyrosine, clone P-TYR-01 (Monoclonal)

Clone no. P-TYR-01

MONOSAN

---

Additional info

-

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES