

Product datasheet MON2058

**MONOSAN**<sup>®</sup>

Mouse anti Hepatocytes, clone EP-1

Clone no. EP-1

MONOSAN

---

Product name	Mouse anti Hepatocytes, clone EP-1
Host	Mouse
Applications	FC, IHC-fr, IF
Species reactivity	human
Conjugate	-
Immunogen	SK-H1A9-2 human hepatocarcinoma cells
Isotype	IgG2b
Clonality	Monoclonal
Clone number	EP-1
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti Hepatocytes, clone EP-1

Clone no. EP-1

MONOSAN

---

**Additional info**

Monoclonal antibody EP-1 stains human liver canaliculi and a subset of hepatocellular carcinomas

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES