

## Mouse anti-Hepatocytes, clone EP-2 (Monoclonal)

Clone no. EP-2

MONOSAN

---

Product name	Mouse anti-Hepatocytes, clone EP-2 (Monoclonal)
Host	Mouse
Applications	FC, IHC-fr, IF
Species reactivity	human
Conjugate	-
Immunogen	HEB-3B human hepatocarcinoma cells
Isotype	IgG2b
Clonality	Monoclonal
Clone number	EP-2
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Hepatocytes, clone EP-2 (Monoclonal)

Clone no. EP-2

MONOSAN

---

**Additional info**

Monoclonal antibody EP-2 stains human hepatocytes and the majority of human hepatocellular carcinomas

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES