Product datasheet MON2071



Mouse anti-CD80, clone MEM-233 (Monoclonal)

Clone no. MEM-233 MONOSAN

Product name Mouse anti-CD80, clone MEM-233 (Monoclonal)

Host Mouse

Applications FC , IP

Species reactivity Human

Conjugate -

Immunogen Extracellular domain of human CD80 fused to human IgG1(Fc)

lsotype lgG1

Clonality Monoclonal

Clone number MEM-233

Size 0.1 mg

Concentration 1 mg/ml

Format -

Storage buffer Phosphate buffered saline (PBS) solution with 15 mM sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON2071



Mouse anti-CD80, clone MEM-233 (Monoclonal)

Clone no. MEM-233 MONOSAN

Additional info

CD80 (B7-1) and CD86 (B7-2) are ligands of T cell critical costimulatory molecule CD28 and of an inhibitory receptor CTLA-4 (CD152). The both B7 molecules are expressed on professional antigen-presenting cells and are essential for T cell activation, the both molecules can also substitute for each other in this process. The question what are the differences in CD80 and CD86 competency has not been fully elucidated yet; there are still conflicts in results about their respective roles in initiation or sustaining of the T cell immune response.

References 1. -

2 -

3.

4. -5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES