

Mouse anti-HLA-G, clone MEM-G/9, Biotin (Monoclonal)

Clone no. MEM-G/9

MONOSAN

Product name	Mouse anti-HLA-G, clone MEM-G/9, Biotin (Monoclonal)
Host	Mouse
Applications	FC , IP, IHC-fr, ICC, ELISA
Species reactivity	Human
Conjugate	Biotin
Immunogen	Recombinant human HLA-G refolded with beta2-microglobulin and peptide.
Isotype	IgG1
Clonality	Monoclonal
Clone number	MEM-G/9
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON2075B

MONOSAN[®]

Mouse anti-HLA-G, clone MEM-G/9, Biotin (Monoclonal)

Clone no. MEM-G/9

MONOSAN

Additional info

-

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES