

Mouse anti-HLA-G, clone MEM-G/1 (Monoclonal)

Clone no. MEM-G/1

MONOSAN

Product name	Mouse anti-HLA-G, clone MEM-G/1 (Monoclonal)
Host	Mouse
Applications	WB, IHC-P (1:60-1:100), IHC-fr
Species reactivity	Human
Conjugate	-
Immunogen	Denatured bacterially expressed recombinant human HLA-G heavy chain.
Isotype	IgG1
Clonality	Monoclonal
Clone number	MEM-G/1
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON2076

MONOSAN[®]

Mouse anti-HLA-G, clone MEM-G/1 (Monoclonal)

Clone no. MEM-G/1

MONOSAN

Additional info

-

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES