

Mouse anti-Placental alkaline phosphatase, clone H7E8 (Monoclonal)

Clone no. H7E8

MONOSAN

Product name	Mouse anti-Placental alkaline phosphatase, clone H7E8 (Monoclonal)
Host	Mouse
Applications	WB, IHC-P, IHC-fr, ELISA
Species reactivity	Human
Conjugate	-
Immunogen	Purified placental alkaline phosphatase.
Isotype	IgG2a
Clonality	Monoclonal
Clone number	H7E8
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON2088

MONOSAN[®]

Mouse anti-Placental alkaline phosphatase, clone H7E8 (Monoclonal)

Clone no. H7E8

MONOSAN

Additional info

-

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES