

Mouse anti-TCR, clone MEM-262 (Monoclonal)

Clone no. MEM-262

MONOSAN

Product name	Mouse anti-TCR, clone MEM-262 (Monoclonal)
Host	Mouse
Applications	FC , IP, WB
Species reactivity	Human
Conjugate	-
Immunogen	Human thymoma cell line HPB-ALL.
Isotype	IgG2a
Clonality	Monoclonal
Clone number	MEM-262
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

