

Mouse anti-Pan-Keratin, clone 80 (Monoclonal)

Clone no. 80

MONOSAN

Product name	Mouse anti-Pan-Keratin, clone 80 (Monoclonal)
Host	Mouse
Applications	IHC-fr (1:20), WB, FC
Species reactivity	human, mouse, rat
Conjugate	-
Immunogen	Callus cytokeratins from human skin
Isotype	IgG1
Clonality	Monoclonal
Clone number	80
Size	1 ml
Concentration	15 ug/ml
Format	-
Storage buffer	Culture medium with 0.01M Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Pan-Keratin, clone 80 (Monoclonal)

Clone no. 80

MONOSAN

Additional info

The antibody does stain all types of keratin containing (that is epithelial) cells in frozen sections of various tissues, with the exception of myo-epithelial cells and the podocytes of the glomeruli.

References

1. Muijen GNP van (1984) Am J Pathol 114, 9
2. Muijen GNP van (1984) Am J. Pathol 116, 363
3. Muijen GNP van (1985) Hum Pathol 16, 590
4. Corver WE et al. (2000) Cytometry 41, 73-80
5. Broek, JCM van et al. (2000) Appl. Immunohis. & Mol. Morph. 8, 316-321

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES