

Mouse anti-Neurofilament, clone 2F11 (Monoclonal)

Clone no. 2F11

MONOSAN

Product name	Mouse anti-Neurofilament, clone 2F11 (Monoclonal)
Host	Mouse
Applications	IHC-P (1:50), IHC-fr (1:50), WB
Species reactivity	human
Conjugate	-
Immunogen	Neurofilament from human brain
Isotype	IgG1
Clonality	Monoclonal
Clone number	2F11
Size	1 ml
Concentration	15 ug/ml
Format	-
Storage buffer	Culture medium with 0.01M Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Neurofilament, clone 2F11 (Monoclonal)

Clone no. 2F11

MONOSAN

Additional info

The antibody does stain neurons in sections of brain and other tissues. On an immunoblot on a crude human brain extract the antibody will only stain the 70 kDa band of neurofilament (ref. 1)

References

1. Muijen GNP van et al. (1984) Am J Pathol. 116, 363-369
2. Klück P et al. (1984) Lancet (March 24th), 652-653
3. Muijen GNP van et al. (1985) Hum Pathol 16, 590-595
4. Banks-Schlegel SP et al. (1985) Cancer Res 406, 339
5. Hoefler H et al. (1985). Virch Arch 406, 339

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES