

Mouse anti-Keratin 5 and 18, clone C-50 (Monoclonal)

Clone no. C-50

MONOSAN

---

Product name	Mouse anti-Keratin 5 and 18, clone C-50 (Monoclonal)
Host	Mouse
Applications	IHC-P, IP, WB, ICC
Species reactivity	Rat, Hamster, Pig, Sheep, Dog, Cow, Human, Mouse
Conjugate	-
Immunogen	Cytoskeleton preparation from HeLa cells.
Isotype	IgG1
Clonality	Monoclonal
Clone number	C-50
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Keratin 5 and 18, clone C-50 (Monoclonal)

Clone no. C-50

MONOSAN

**Additional info**

Cytokeratins are a subfamily of intermediate filaments and characterized by remarkable biochemical diversity. They are represented in epithelial tissues by at least 20 different polypeptides, molecular weight between 40 kDa and 68 kDa. The individual cytokeratin polypeptides are designated 1 to 20 and divided into the type I (acidic cytokeratins 9-20) and type II (basic to neutral cytokeratins 1-8) families.

**References**

1. -
2. -
3. -
4. -
5. -

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**