

Mouse anti-Histone H1, clone AE-4 (Monoclonal)

Clone no. AE-4

MONOSAN

Product name	Mouse anti-Histone H1, clone AE-4 (Monoclonal)
Host	Mouse
Applications	FC, IHC-fr, IHC-P, IF
Species reactivity	human, mouse, rat
Conjugate	-
Immunogen	human leukemia biopsy cells
Isotype	IgG1
Clonality	Monoclonal
Clone number	AE-4
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Histone H1, clone AE-4 (Monoclonal)

Clone no. AE-4

MONOSAN

Additional info

AE-4 recognizes histone H1 in the nuclei of animal and plant cells. It produces a speckled nuclear staining pattern in normal and malignant cells.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES