

Mouse anti-CA-19.9, clone 121SLE (Monoclonal)

Clone no. 121SLE

MONOSAN

---

Product name	Mouse anti-CA-19.9, clone 121SLE (Monoclonal)
Host	Mouse
Applications	IHC-P, ELISA
Species reactivity	Human
Conjugate	-
Immunogen	Immunoprecipitate obtained after immunodiffusion of MAb 19-9 and mucins isolated from an ovarian cyst of a 0Le(a+b-) patient.
Isotype	IgM
Clonality	Monoclonal
Clone number	121SLE
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Tris buffered saline (TBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON3060

**MONOSAN**<sup>®</sup>

Mouse anti-CA-19.9, clone 121SLE (Monoclonal)

Clone no. 121SLE

MONOSAN

---

Additional info

-

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

---