Product datasheet MON3079



Mouse anti-Lck-interacting molecule, clone LIME-06 (Monoclonal)

Clone no. LIME-06 MONOSAN

Product name Mouse anti-Lck-interacting molecule, clone LIME-06 (Monoclonal)

Host Mouse

Applications IP, FC

Species reactivity Human

Conjugate -

Immunogen Bacterially expressed intracellular fragment corresponding to aa 141-295 of

human LIME.

lsotype lgG1

Clonality Monoclonal

Clone number LIME-06

Size 0.1 mg

Concentration 1 mg/ml

Format -

Storage buffer Phosphate buffered saline (PBS) solution with 15 mM sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON3079



Mouse anti-Lck-interacting molecule, clone LIME-06 (Monoclonal)

Clone no. LIME-06 MONOSAN

Additional info

LIME (Lck-interacting molecule) is a 30 kDa double-palmitoylated protein with unusually basic cytoplasmic domain, expressed by T cells. After ligation of CD4 or CD8 T cell coreceptors, LIME is phosphorylated by Src-family kinases and associates with Lck and Fyn kinases and with their negative regulator Csk. Interestingly, Csk-mediated phosphorylation of C-terminal negative-regulatory tyrosine of LIME-associated Lck can result in increase of enzymatic activity compared with the total pool of Lck, thus, LIME serves as a positive regulator of TCR-dependent T cell signaling. However, under some circumstances, LIME may mediate inhibitory signals.

References 1.

2 -

3. -

4. -

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES