

Mouse anti-Cytokeratin 19, clone A53-B/A2.26 (Monoclonal)

Clone no. A53-B/A2.26

MONOSAN

---

Product name	Mouse anti-Cytokeratin 19, clone A53-B/A2.26 (Monoclonal)
Host	Mouse
Applications	IHC-P (1:100-1:500)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	A53-B/A2.26
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Cytokeratin 19, clone A53-B/A2.26 (Monoclonal)

Clone no. A53-B/A2.26

MONOSAN

**Additional info**

Anti-Cytokeratin 19 reacts with a wide variety of epithelium and epithelial malignancies including Adenocarcinomas of the Colon, stomach, pancreas, biliary tract, liver and breast. Perhaps the most useful application is the identification of Thyroid carcinoma of the Papillary type, although Follicular Carcinoma is also labeled by this antibody approximately 50-60% of the time.

**References**

1. Cerilli LA, et al. Am J Clin Pathol 2002;118:186-193
2. Jain R, et al. Appl Immunohistochem Mol Morphol. 2010; 18:9-15
3. Rosai J. Tumori. 2003; 89:517-9
4. de Matos LL, et al. Diagn Pathol. 2012; 7:97
5. Zubair W, et al. Hum Pathol 30:1166-1171

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**