

Mouse anti-Human PSAP, clone PASE/4LJ (Monoclonal)

Clone no. PASE/4LJ

MONOSAN

Product name	Mouse anti-Human PSAP, clone PASE/4LJ (Monoclonal)
Host	Mouse
Applications	IHC-P (1:25-1:100)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	PASE/4LJ
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human PSAP, clone PASE/4LJ (Monoclonal)

Clone no. PASE/4LJ

MONOSAN

Additional info

The antibody reacts with prostatic acid phosphatase in the glandular epithelium of normal and hyperplastic prostate, and adenocarcinoma of the prostate. Anti-PSAP is useful in identifying prostatic origin of tumors in the metastatic setting. PSAP complements other immunohistochemical markers in the correct clinical context.

References

1. Ansari, MA, et al. Am J Clin Path 1981;76:94-98
2. Kimura, N, et al. Virchows Arch A 1986;4:247-251
3. Kidwai N et al. Breast Cancer Res. 2004;6(1):R18-23
4. Genega EM et al. Mod Pathol. 2000 Nov;13(11):1186-91
5. Gatalica Z et al. Appl Immunohistochem Mol Morphol. 2000 Jun;8(2):158-61

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES