

Mouse anti-Human Renal Cell Carcinoma, clone PN-15 (Monoclonal)

Clone no. PN-15

MONOSAN

Product name	Mouse anti-Human Renal Cell Carcinoma, clone PN-15 (Monoclonal)
Host	Mouse
Applications	IHC-P (1:50-1:200)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1-k
Clonality	Monoclonal
Clone number	PN-15
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human Renal Cell Carcinoma, clone PN-15 (Monoclonal)

Clone no. PN-15

MONOSAN

Additional info

Anti-renal cell carcinoma (RCC) recognizes a 200 kD glycoprotein localized in the brush border of the proximal renal tubule. This antibody immunoreacts with approximately most primary renal cell carcinomas and can aid in the diagnosis when renal cell carcinoma enters the differential diagnosis. Other tumors that may react with this antibody are parathyroid adenoma, an occasional breast carcinoma. Nephroblastoma, oncocytoma, mesoblastic nephroma, transitional cell carcinoma, and angiomyolipoma are not labeled with this antibody.

References

1. Dabbs, D. 4th Edition. Elsevier Saunders. 2014; p234
2. Bakshi N, et al. Appl Immunohistochem Mol Morphol. 2007; 15:310-5
3. McGregor DK, et al. Am J Surg Pathol. 2001; 25:1485-92
4. Avery, AK et al. Am J Surg Pathol 24(2): 203-210, 2000
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES