

Mouse anti-Human Spectrin, clone RBC2/3D5 (Monoclonal)

Clone no. RBC2/3D5

MONOSAN

Product name	Mouse anti-Human Spectrin, clone RBC2/3D5 (Monoclonal)
Host	Mouse
Applications	IHC-P (1:10-1:100)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG2b-k
Clonality	Monoclonal
Clone number	RBC2/3D5
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human Spectrin, clone RBC2/3D5 (Monoclonal)

Clone no. RBC2/3D5

MONOSAN

Additional info

Spectrin is a cytoskeletal protein which is found in muscles, red blood cells and red cell precursors. Anti-Spectrin antibody is useful in the identification of blood dyscrasias and muscle disorders.

References

1. Sadahira, Y et al. J Clin Pathol. 1999 Dec; 52(12): 919-21
2. Nehls, V et al. Am J Pathol. 1993 May; 142(5): 1565-73
3. Muller M et al. J Vet Med A Physiol Pathol Clin Med. 2001 Feb;48(1):51-7
4. Terada N et al. J Anat. 1997 Apr;190(Pt 3):397-404
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES