

Rabbit anti-Human Epstein-Barr Virus, clone MRQ-47 (Monoclonal)

Clone no. MRQ-47

MONOSAN

Product name	Rabbit anti-Human Epstein-Barr Virus, clone MRQ-47 (Monoclonal)
Host	Rabbit
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	MRQ-47
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rabbit anti-Human Epstein-Barr Virus, clone MRQ-47 (Monoclonal)

Clone no. MRQ-47

MONOSAN

Additional info

Analyte Specific Reagent: Analytical and performance characteristics are not determined.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES