Product datasheet MON3397



Mouse anti-MUC4, clone 8G7 (Monoclonal)

Clone no. 8G7 MONOSAN

Product name Mouse anti-MUC4, clone 8G7 (Monoclonal)

Host Mouse

Applications IHC-P (1:10-1:50)

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype IgG1-k

Clonality Monoclonal

Clone number 8G7

Size 1 ml

Concentration n/a

Format -

Storage buffer Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON3397



Mouse anti-MUC4, clone 8G7 (Monoclonal)

Clone no. 8G7 MONOSAN

Additional info

MUC4 (mucin 4) is a large membrane-anchored glycoprotein of the mucin family that is expressed by epithelial cells in various normal tissues including lung, bronchus, stomach, colon, and cervix. MUC4 is generally not detected in normal pancreas, but is expressed in the vast majority of pancreatic neoplasms, such as pancreatic ductal adenocarcinoma. Additionally, expression in various carcinomas has been reported, including gastric adenocarcinoma, colon adenocarcinoma, and lung adenocarcinoma.

References 1. Moniaux N, et al. J Histochem Cytochem. 2004; 52:253-61

2 Moniaux N, et al. Br J Cancer. 2004;91:1633-8

3. Llinares K, et al. Mod Pathol. 2004;17:150-7

4. -5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES