Product datasheet

MON3402



MONOSAN

Mouse anti-DC-SIGN/CD209, clone DCN47.5 (Monoclonal) Clone no. DCN47.5

Product name	Mouse anti-DC-SIGN/CD209, clone DCN47.5 (Monoclonal)
Host	Mouse
Applications	FC,FUNC
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietery to MONOSAN and/or its suppliers
lsotype	lgG1
Clonality	Monoclonal
Clone number	DCN47.5
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON3402

MONOSAN

Mouse anti-DC-SIGN/CD209, clone DCN47.5 (Monoclonal) Clone no. DCN47.5

MONOSAN

Additional info

The monoclonal antibody DCN47.5 reacts with the C-type lectin, DC-SIGN (CD209), exclusively expressed on human dendritic cells (DC). DC are specialized antigen presenting cells and bridge the innate and the adaptive immune system. They provide high levels of costimulation necessary for activation of both naà ve and antigen-experienced T-cells. Immature DC are capable to migrate to inflammatory sites, capture antigen and process it internally to form MHC-peptide complexes. Following antigen uptake, DC undergo maturation and migrate to lymphoid organs where they can present MHC-peptide complexes to resting T-cells and drive T-cell proliferation. During differentiation and maturation of DC, several phenotypic surface markers are expressed: CD1a, CD4, CD11, CD40, CD86, and HLA-DR. Immature DC predominantly express CCR5 which enables DC to migrate to inflammatory sites, whereas mature DC express high levels of CXCR4, a receptor that facilitates migration to lymphoid organs.</br> DC-specific, ICAM-3 grabbing, nonintegrin (DC-SIGN), a 44 kDa C-type lectin that binds to the HIV-1 envelope glycoprotein gp120, ICAM-3 on T-cells and ICAM-2 on endothelial cells. HIV virions are able to infect cells expressing CD4 and the chemokine co-receptors CCR5 or CXCR4 and can attach to DC-SIGN to extend virion lifespan. Therefore, DC are candidates for HIV-1 infection. DC-SIGN-ICAM-3 binding is integrin-independent but calciumdependent and antibodies reactive against DC-SIGN can be used to study DC-SIGN-ICAM3 binding.</br> reacts with the C-type lectin DC-SIGN (CD209) expressed on human dendritic cells and inhibits binding of DC-SIGN to ICAM-2 on endothelial cells.

References	1.	Leeuwenberg; JFM et al. J Immunol 1994; 152: 4036
	2	Leeuwenberg, JFM et al J Immnol 1994, 152: 5070
	3.	Marchetti; L et al. J Biol Chem 2004; 279: 32869
	4.	-
	5	_

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com

FOR-044 22-12-2021