

Mouse anti-Actin, clone 647 (Monoclonal)

Clone no. 647

MONOSAN

Product name	Mouse anti-Actin, clone 647 (Monoclonal)
Host	Mouse
Applications	IHC-fr, IHC-P (1:5-1:10), IF, WB
Species reactivity	Chicken, Human, Mouse, Rat
Conjugate	-
Immunogen	Antigen isolated from chicken skeletal muscle cells
Isotype	IgM
Clonality	Monoclonal
Clone number	647
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Contains 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Actin, clone 647 (Monoclonal)

Clone no.

647

MONOSAN

Additional info

Mab 647 is useful for studying the intracellular distribution and structure of actin in the cytoskeletal system. In immunoblotting one single band is reactive corresponding to actin.

Positive control: Muscle or myofibrils.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES