

## Mouse anti-Myelin, clone 2B5 (Monoclonal)

Clone no. 2B5

MONOSAN

---

Product name	Mouse anti-Myelin, clone 2B5 (Monoclonal)
Host	Mouse
Applications	IHC-fr, IHC-P (1:10), WB
Species reactivity	Human
Conjugate	-
Immunogen	Human myelin isolated from peripheral nervous system
Isotype	IgM
Clonality	Monoclonal
Clone number	2B5
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Contains 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

## Mouse anti-Myelin, clone 2B5 (Monoclonal)

Clone no. 2B5

MONOSAN

**Additional info**

In paraffin sections antibody pm 43 recognizes only myelin in the peripheral nervous system (PNS). In frozen sections it also reacts with myelin in the central nervous system (CNS). Is useful for studying myelination and demyelination processes in the PNS, and remyelination processes in the CNS as can be observed after trauma in multiple sclerosis. Reacts with a 43 kD protein in immunoblotting.

Positive control: Peripheral nerve.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES