Product datasheet MON4021



Mouse anti-SHP2, clone SIT-01 (Monoclonal)

Clone no. SIT-01 MONOSAN

Product name Mouse anti-SHP2, clone SIT-01 (Monoclonal)

Host Mouse

Applications FC, IP, WB

Species reactivity Human

Conjugate -

Immunogen Bacterially produced recombinant intracellular fragment of human SIT.

lsotype lgG1

Clonality Monoclonal

Clone number SIT-01

Size 0.1 mg

Concentration 1 mg/ml

Format -

Storage buffer Phosphate buffered saline (PBS) solution with 15 mM sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON4021



Mouse anti-SHP2, clone SIT-01 (Monoclonal)

Clone no. SIT-01 MONOSAN

Additional info

SIT (SHP2-interacting transmembrane adaptor protein) is expressed exclusively in lymphoid organs and acts either as a positive or as a negative regulatory element in T cell activation and in T cell development. Binding to Grb2 plays a pivotal role in signal transduction. Hubener et al. (2001) determined that the SIT gene contains 5 exons and spans 1.8 kb of genomic DNA. The SIT promoter demonstrated strong transcriptional activity and potential binding sites for both ubiquitous and lymphoid-specific transcription factors.

References 1. -

2 -

3. -

4. -

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES