

Mouse anti-TLR2, clone TL2.3, Biotin (Monoclonal)

Clone no. TL2.3

MONOSAN

---

Product name	Mouse anti-TLR2, clone TL2.3, Biotin (Monoclonal)
Host	Mouse
Applications	-
Species reactivity	-
Conjugate	Biotin
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Monoclonal
Clone number	TL2.3
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-TLR2, clone TL2.3, Biotin (Monoclonal)

Clone no. TL2.3

MONOSAN

Additional info

-

## References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES