Product datasheet MON4032B



## Rat anti Mouse-TLR4/MD2, clone MTS510, Biotin (Monoclonal)

Clone no. MTS510 MONOSAN

Product name Rat anti Mouse-TLR4/MD2, clone MTS510, Biotin (Monoclonal)

**Host** Rat

Applications -

Species reactivity -

**Conjugate** Biotin

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype -

**Clonality** Monoclonal

Clone number MTS510

Size 1 ml

**Concentration** n/a

Format -

Storage buffer PBS with 0.1% BSA and 0.02% sodium azide

Storage until expiry date 2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON4032B



Rat anti Mouse-TLR4/MD2, clone MTS510, Biotin (Monoclonal)

Clone no. MTS510 MONOSAN

Additional info

References 1. -

2 -

3.

4. -

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES