

Mouse anti-HIV p24 protein, clone ND1 (Monoclonal)

Clone no. ND1

MONOSAN

Product name	Mouse anti-HIV p24 protein, clone ND1 (Monoclonal)
Host	Mouse
Applications	WB, ICC, ELISA
Species reactivity	Human Immunodeficiency Virus
Conjugate	-
Immunogen	HIV-1 p24 C-terminal peptide.
Isotype	IgG1
Clonality	Monoclonal
Clone number	ND1
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-HIV p24 protein, clone ND1 (Monoclonal)

Clone no. ND1

MONOSAN

Additional info

-

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES