

Mouse anti-CD82, clone B-L2 (Monoclonal)

Clone no. B-L2

MONOSAN

Product name	Mouse anti-CD82, clone B-L2 (Monoclonal)
Host	Mouse
Applications	FC,IHC-fr,IF
Species reactivity	human
Conjugate	Purified
Immunogen	CD82 transfected murine mastocytoma cell line P815.
Isotype	IgG1
Clonality	Monoclonal
Clone number	B-L2
Size	2 ml
Concentration	n/a
Format	-
Storage buffer	PBS with 0.09% sodium azide and 1% BSA
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD82, clone B-L2 (Monoclonal)

Clone no. B-L2

MONOSAN

Additional info

Mouse anti Human CD82 antibody, clone B-L2 recognizes human CD82, also known as C33 antigen, Metastasis suppressor Kangai-1, Suppressor of tumorigenicity 6 protein or Tetraspanin-27. CD82 is a 267 amino acid ~50-53 kDa multiple pass transmembrane glycoprotein and member of the tetraspanin family expressed by T cells, NK cells, monocytes, granulocytes and platelets.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES