

Mouse anti-CD103, clone LF61, R-PE (Monoclonal)

Clone no. LF61

MONOSAN

Product name	Mouse anti-CD103, clone LF61, R-PE (Monoclonal)
Host	Mouse
Applications	FC
Species reactivity	human
Conjugate	RPE
Immunogen	Hairy cell leukemia cells
Isotype	IgG1
Clonality	Monoclonal
Clone number	LF61
Size	100 TESTS
Concentration	n/a
Format	-
Storage buffer	PBS with 0.09% sodium azide and 1% BSA, 5% sucrose
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD103, clone LF61, R-PE (Monoclonal)

Clone no. LF61

MONOSAN

Additional info

Mouse anti Human CD103 antibody, clone LF61 recognizes the human CD103 cell surface antigen, a glycoprotein expressed by approximately 1% of peripheral blood lymphocytes, activated T lymphocytes and by hairy cell leukemia cells. The antigen is also expressed by intraepithelial lymphocytes. It has recently been shown to be identical to the alpha E integrin.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES