

Mouse anti-CD89, clone MIP8a, R-PE (Monoclonal)

Clone no. MIP8a

MONOSAN

---

Product name	Mouse anti-CD89, clone MIP8a, R-PE (Monoclonal)
Host	Mouse
Applications	FC
Species reactivity	human
Conjugate	RPE
Immunogen	Recombinant soluble human Fc alpha R.
Isotype	IgG1
Clonality	Monoclonal
Clone number	MIP8a
Size	100 TESTS
Concentration	n/a
Format	-
Storage buffer	PBS with 0.09% sodium azide and 1% BSA
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD89, clone MIP8a, R-PE (Monoclonal)

Clone no. MIP8a

MONOSAN

**Additional info**

Mouse anti Human CD89 antibody, clone MIP8a recognizes the human CD89 cell surface antigen, a 50-75 kDa cell surface glycoprotein that is also known as the IgA receptor (Fc alpha R). CD89 is expressed by peripheral blood monocytes and neutrophils. MIP8a blocks binding of IgA to the Fc alpha R, and also inhibits neutrophil phagocytosis of IgA complexes.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES