Product datasheet MON4073



Mouse anti-CD200, clone OX-104 (Monoclonal)

Clone no. OX-104 MONOSAN

Product name Mouse anti-CD200, clone OX-104 (Monoclonal)

Host Mouse

Applications FC,IHC-fr

Species reactivity human

Conjugate Purified

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

lsotype lgG1

Clonality Monoclonal

Clone number OX-104

Size 0.1 mg

Concentration 1.0mg/ml

Format -

Storage buffer PBS with 0.09% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON4073



Mouse anti-CD200, clone OX-104 (Monoclonal)

Clone no. OX-104 MONOSAN

Additional info

Mouse anti Human CD200 antibody, clone OX-104 recognizes the human CD200 cell surface antigen, also known as OX2.CD200 is expressed by a subset of B lymphocytes, some endothelial cells and by neurons. Studies have suggested that the CD200-CD200 ligand system is of importance in the control of macrophage and granulocyte activation.

References 1. -

2 -

3. -

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES