Product datasheet

MON4078



MONOSAN

Mouse anti-CD173, clone BRIC231 (Monoclonal) Clone no. BRIC231

Product name	Mouse anti-CD173, clone BRIC231 (Monoclonal)	
Host	Mouse	
Applications	IF,FC,FUNC	
Species reactivity	human,pig	
Conjugate	Purified	
Immunogen	Human erythroleukemic cell line (HEL) established from a 30 year old patient with relapsed erythroleukemia following treatment for Hodgkin lymphoma.	
lsotype	lgG1	
Clonality	Monoclonal	
Clone number	BRIC231	
Size	0.2 mg	
Concentration	1.0mg/ml	
Format	-	
Storage buffer	TRIS buffered glycine	
Storage until expiry date	2-8°C	

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON4078



Mouse anti-CD173, clone BRIC231 (Monoclonal) Clone no. BRIC231

MONOSAN

Additional info

Mouse anti Human CD173 antibody, clone BRIC231 recognizes human type 2 H blood group antigen, also known as CD173. Active H substances in man, are expressed by many cells and tissues and also by erythrocytes.

References	1.	-
	2	-
	3.	-
	4.	-
	5.	-

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com