

Mouse anti-CD88, clone P12/1 (Monoclonal)

Clone no. P12/1

MONOSAN

---

Product name	Mouse anti-CD88, clone P12/1 (Monoclonal)
Host	Mouse
Applications	WB,IHC-P,FC,IHC-fr
Species reactivity	human,rhesus monkey
Conjugate	Purified
Immunogen	C5aR - peptide: Met1- Asn31
Isotype	IgG2a
Clonality	Monoclonal
Clone number	P12/1
Size	0.2 mg
Concentration	1 mg/ml
Format	-
Storage buffer	PBS with 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD88, clone P12/1 (Monoclonal)

Clone no. P12/1

MONOSAN

**Additional info**

Mouse anti Human CD88 antibody, clone P12/1 recognizes the C5a receptor (C5aR) also known as CD88 or C5a anaphylatoxin chemotactic receptor 1. CD88 is predominantly expressed on cells of the myeloid lineage. When C5aR is preincubated with C5a, Mouse anti Human CD88 antibody, clone P12/1 does not bind to the receptor, as the binding site of P12/1 is located in the C5a binding region (Werfel et al. 1996 and Weinman et al. 2003)

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES