Product datasheet

Product name

Host

MON4091



MONOSAN

Rat anti-CD195, clone HEK/1/85a (Monoclonal) HEK/1/85a Clone no.

Rat anti-CD195, clone HEK/1/85a (Monoclonal)
Rat

Applications	FC
Species reactivity	human
Conjugate	Purified
Immunogen	CHO cells transfected with human CCR5
lsotype	lgG2a
Clonality	Monoclonal
Clone number	HEK/1/85a
Size	0.1 mg
Concentration	1.0 mg/ml
Format	-
Storage buffer	PBS with 0.09% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON4091



Rat anti-CD195, clone HEK/1/85a (Monoclonal) Clone no. HEK/1/85a

MONOSAN

Additional info

Rat anti Human CD195 antibody, clone HEK/1/85a recognizes the human CD195 cell surface antigen, a 45 kDa glycoprotein also known as CCR5.CD195 acts as a receptor for a number of chemokines including RANTES and eotaxin and serves as a co-receptor for the entry of HIV into cells. CD195 is expressed by a subset of T lymphocytes and by monocytes.

References	1.	-
	2	-
	3.	-
	4.	-
	5.	-

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com