Product datasheet MON4091R



## Rat anti-CD195, clone HEK/1/85a, R-PE (Monoclonal)

Clone no. HEK/1/85a MONOSAN

Product name Rat anti-CD195, clone HEK/1/85a, R-PE (Monoclonal)

**Host** Rat

**Applications** FC

Species reactivity human

**Conjugate** RPE

Immunogen CHO cells transfected with human CCR5

**Isotype** IgG2a

**Clonality** Monoclonal

Clone number HEK/1/85a

Size 100 TESTS

**Concentration** n/a

Format -

Storage buffer PBS with 0.09% sodium azide and 1% BSA

Storage until expiry date 2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON4091R



Rat anti-CD195, clone HEK/1/85a, R-PE (Monoclonal)

Clone no. HEK/1/85a MONOSAN

## Additional info

Rat anti Human CD195 antibody, clone HEK/1/85a recognizes the human CD195 cell surface antigen, a 45 kDa glycoprotein also known as CCR5.CD195 acts as a receptor for a number of chemokines including RANTES and eotaxin and serves as a co-receptor for the entry of HIV into cells. CD195 is expressed by a subset of T lymphocytes and by monocytes.

References 1. -

2 -

3.

4. -

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES