Product datasheet MON4098



Mouse anti-CD226, clone DX11 (Monoclonal)

Clone no. DX11 MONOSAN

Product name Mouse anti-CD226, clone DX11 (Monoclonal)

Host Mouse

Applications IP,FC

Species reactivity human

Conjugate Purified

Immunogen Human cytotoxic T lymphocyte clone

lsotype lgG1

Clonality Monoclonal

Clone number DX11

Size 0.1 mg

Concentration 1.0 mg/ml

Format -

Storage buffer PBS with 0.09% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON4098



Mouse anti-CD226, clone DX11 (Monoclonal)

Clone no. DX11 MONOSAN

Additional info

Mouse anti Human CD226 antibody, clone DX11 recognizes human CD226, a ~65 kDa glycoprotein, also known as DNAM1 (DNAX accessory molecule-1). CD226 is broadly expressed on T-cells, NK cells, platelets, monocytes and a subset of B cells. CD226 is also expressed by a subset of CD3 positive thymocytes. Mouse anti Human CD226 antibody, clone DX11 is reported to inhibit T- and NK cell mediated cytotoxicity against tumor cell targets and to block TNF alpha and IFN gamma secretion by alloantigen-specific T-cells (Kojima et al. 2003).

References 1.

2 -

3. -

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES