Product datasheet MON4104S



## Mouse anti-Giardia Lamblia, clone 6231

Clone no. 6231 MONOSAN

Product name Mouse anti-Giardia Lamblia, clone 6231

**Host** Mouse

**Applications** WB,IF

Species reactivity protozoan

Conjugate Purified

Immunogen 65 kDa <em>Giardia</em> antigen.

lsotype lgG1

**Clonality** Monoclonal

Clone number 6231

Size 0.1 mg

Concentration 0.1mg/ml

Format -

Storage buffer PBS with 0.09% sodium azide

Storage until expiry date 2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON4104S



Mouse anti-Giardia Lamblia, clone 6231

Clone no. 6231 MONOSAN

## Additional info

Mouse anti Giardia lamblia antibody, clone 6231 recognizes Giardia lamblia, a flagellated protozoan parasite which infects the small intestines of humans and a variety of other mammalian hosts. Infection by G. lamblia results in giardiasis; a type of gastroenteritis that causes severe diarrhea and abdominal cramps. It has been suggested that chronic giardiasis can result in long-term growth retardation (Simsek. et al. 2004).

References 1. -

2 -

3.

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES