

Mouse anti-Human kappa light chain, clone 2B7 (Monoclonal)

Clone no. 2B7

MONOSAN

---

Product name	Mouse anti-Human kappa light chain, clone 2B7 (Monoclonal)
Host	Mouse
Applications	IHC-Fr (1:30), WB
Species reactivity	human
Conjugate	-
Immunogen	kappa Bence Jones protein
Isotype	IgG2b
Clonality	Monoclonal
Clone number	2B7
Size	1 ml
Concentration	10 ug/ml
Format	-
Storage buffer	Culture medium with 0.01M Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human kappa light chain, clone 2B7 (Monoclonal)

Clone no. 2B7

MONOSAN

---

**Additional info**

The antibody reacts with free kappa chains, but not with kappa chains attached to heavy chains. No reaction is observed with Ig kappa on membranes of lymphocytes. No crossreactivity with lambda chains. In Ouchterlony immunodiffusion no reaction with Bence Jones kappa is observed.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES