

Mouse anti-Beta-Galactosidase from E. coli, clone BG-02 (Monoclonal)

Clone no. BG-02

MONOSAN

Product name	Mouse anti-Beta-Galactosidase from E. coli, clone BG-02 (Monoclonal)
Host	Mouse
Applications	WB, ICC
Species reactivity	Escherichia coli, Tagged fusion proteins in all species
Conjugate	-
Immunogen	beta-Galactosidase purified from <i>E. coli</i>.
Isotype	IgG1
Clonality	Monoclonal
Clone number	BG-02
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Beta-Galactosidase from E. coli, clone BG-02 (Monoclonal)

Clone no. BG-02

MONOSAN

Additional info

-

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES