

Mouse anti-Interleukin 6, clone 5E1 (Monoclonal)

Clone no. 5E1

MONOSAN

Product name	Mouse anti-Interleukin 6, clone 5E1 (Monoclonal)
Host	Mouse
Applications	FUNC,ELISA,IP,WB
Species reactivity	human, cynomolgus monkey, rhesus monkey
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	5E1
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Interleukin 6, clone 5E1 (Monoclonal)

Clone no. 5E1

MONOSAN

Additional info

The antibody reacts with human native and recombinant IL-6 as assessed by ELISA. The antibody inhibits the biological activity of human native and recombinant IL-6 as determined with the B9 cell bioassay. The antibody cross reacts with rhesus and cynomolgus natural IL6.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES