Product datasheet MON5042



Mouse anti-Biotin, clone BIO-8 (Monoclonal)

Clone no. BIO-8 MONOSAN

Product name Mouse anti-Biotin, clone BIO-8 (Monoclonal)

Host Mouse

Applications IHC, ELISA, FC

Species reactivity tag

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

lsotype lgG3

Clonality Monoclonal

Clone number BIO-8

Size 1 ml

Concentration 100 ug/ml

Format -

Storage buffer Tissue culture supernatant with sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON5042



Mouse anti-Biotin, clone BIO-8 (Monoclonal)

Clone no. BIO-8 MONOSAN

Additional info

The antibody recognizes an epitope that is exposed by protein-biotin and nucleic acid-biotin conjugates. The antibody can be used in molecular biological tests and immunohistology. The antibody permits the formation of antibody-biotin-streptavidin-enzyme complexes thus enhancing the sensitivity of the detection system.

References 1. -

2 -

3. -

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES