

Mouse anti-Streptavidin, clone STREP-10 (Monoclonal)

Clone no. STREP-10

MONOSAN

Product name	Mouse anti-Streptavidin, clone STREP-10 (Monoclonal)
Host	Mouse
Applications	IHC (1:20), ELISA
Species reactivity	tag
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	STREP-10
Size	1 ml
Concentration	200 ug/ml
Format	-
Storage buffer	Supernatant with 1% BSA and 0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Streptavidin, clone STREP-10 (Monoclonal)

Clone no. STREP-10

MONOSAN

Additional info

The antibody is specific for streptavidin. It lacks crossreactivity with avidin.
The antibody does not interfere with biotin binding.

References

1. Scheffer GL et al.
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES