

Rat anti-Mouse IgA, clone LO-MA-7 (Monoclonal)

Clone no. LO-MA-7

MONOSAN

Product name	Rat anti-Mouse IgA, clone LO-MA-7 (Monoclonal)
Host	Rat
Applications	IHC-fr, ELISA
Species reactivity	mouse
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgMk
Clonality	Monoclonal
Clone number	LO-MA-7
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	PBS with 0.1% sodium azide
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-Mouse IgA, clone LO-MA-7 (Monoclonal)

Clone no. LO-MA-7

MONOSAN

Additional info

Mouse alpha heavy chain of immunoglobulin (determined by immunodot).
Avidity on IgA: $4.8 \times 10^9 \text{ M}^{-1}$. Useful as second antibody in immunoassays, such as ELISA's (good binding on plastic) and immunohistochemistry on frozen sections. Available in unconjugated MON5050, HRP conjugated MON5050P, Biotin conjugated MON5050B

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES