

Rat anti-Mouse IgA, clone LO-MA-7, HRP (Monoclonal)

Clone no. LO-MA-7

MONOSAN

Product name	Rat anti-Mouse IgA, clone LO-MA-7, HRP (Monoclonal)
Host	Rat
Applications	IHC-fr, ELISA
Species reactivity	mouse
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgMk
Clonality	Monoclonal
Clone number	LO-MA-7
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	PBS with 0.1% sodium azide 50% glycerol
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-Mouse IgA, clone LO-MA-7, HRP (Monoclonal)

Clone no. LO-MA-7

MONOSAN

Additional info

Mouse alpha heavy chain of immunoglobulin (determined by immunodot).
Avidity on IgA: $2.9 \times 10^9 \text{ M}^{-1}$. Available as unconjugated MON5050, Biotin
conjugated MON5050B, Peroxidase conjugated MON5050P.

References

1. Bazin et al. Pergamon Press Oxford and NY 1982; 29:615-618
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES