

Rat anti-Mouse IgE, clone LO-ME-3, HRP (Monoclonal)

Clone no. LO-ME-3

MONOSAN

---

Product name	Rat anti-Mouse IgE, clone LO-ME-3, HRP (Monoclonal)
Host	Rat
Applications	IHC-fr, ELISA
Species reactivity	mouse
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1k
Clonality	Monoclonal
Clone number	LO-ME-3
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	PBS with 0.1% sodium azide
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-Mouse IgE, clone LO-ME-3, HRP (Monoclonal)

Clone no. LO-ME-3

MONOSAN

---

**Additional info**

Mouse epsilon heavy chain of immunoglobulin (determined by immunodot). Useful as second antibody in immunoassays, such as ELISA's (good binding on plastic) and immunohistochemistry on frozen sections. Available in unconjugated MON5051, FITC conjugated MON5051F, HRP conjugated MON5050P and Biotin conjugated MON5051B.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES