

Rat anti-Mouse IgG1, clone LO-MG1-2, HRP (Monoclonal)

Clone no. LO-MG1-2

MONOSAN

Product name	Rat anti-Mouse IgG1, clone LO-MG1-2, HRP (Monoclonal)
Host	Rat
Applications	IHC-fr, ELISA
Species reactivity	mouse
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgGk
Clonality	Monoclonal
Clone number	LO-MG1-2
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	PBS with 0.1% sodium azide
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-Mouse IgG1, clone LO-MG1-2, HRP (Monoclonal)

Clone no. LO-MG1-2

MONOSAN

Additional info

Mouse gamma 1 heavy chain of immunoglobulin (determined by immunodot). No crossreactivity with other animals. Available as unconjugated MON5053, FITC conjugated MON 5053F, HRP conjugated Cat.no. MON 5053P, Biotin conjugated MON 5053B.

References

1. Bazin, et al., J.Immunology Meth. 1984;71, 9-16
2. Querinjean, et al., Eur.J.Biochem. 1972: 31:354-359
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES