

Rat anti-Mouse IgG2a, clone LO-MG2a-2, Biotin (Monoclonal)

Clone no. LO-MG2a-2

MONOSAN

Product name	Rat anti-Mouse IgG2a, clone LO-MG2a-2, Biotin (Monoclonal)
Host	Rat
Applications	IHC-fr, ELISA
Species reactivity	mouse
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgGk
Clonality	Monoclonal
Clone number	LO-MG2a-2
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	PBS with 0.1% sodium azide
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-Mouse IgG2a, clone LO-MG2a-2, Biotin (Monoclonal)

Clone no. LO-MG2a-2

MONOSAN

Additional info

Mouse Gamma 2a Heavy Chain of Immunoglobulin (determined by immunodot). Avidity on IgG2a: 7×10^9 M⁻¹. No crossreactivity with other animals. Available as unconjugated MON5054, HRP conjugated MON 5054P, Biotin conjugated MON 5054B

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES