

Rat anti-Mouse IgG3, clone LO-MG3-7, HRP (Monoclonal)

Clone no. LO-MG3-7

MONOSAN

---

Product name	Rat anti-Mouse IgG3, clone LO-MG3-7, HRP (Monoclonal)
Host	Rat
Applications	IHC-fr, ELISA
Species reactivity	mouse, swine, horse, dog
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgMk
Clonality	Monoclonal
Clone number	LO-MG3-7
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	PBS with 0.1% sodium azide
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-Mouse IgG3, clone LO-MG3-7, HRP (Monoclonal)

Clone no. LO-MG3-7

MONOSAN

---

**Additional info**

Mouse Gamma 3 Heavy Chain of Immunoglobulin (determined by immunodot). Avidity on IgG3:  $2 \times 10^{10}$  M<sup>-1</sup>. Crossreactivity with: Horse (+), swine (++) and dog (+/-). Available in unconjugated MON5056, FITC conjugated MON5056F, HRP conjugated MON5056H, Biotin conjugated MON5056B.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES